



UNILIFT PRO

Next Generation - Ceiling Lifts

14:02 Uhr
Battery 

76 kg

Features and Options

- > Integrated scale
- > Steel internal parts, no plastic gears
- > Lifting capacity 600 up to 1,100 lbs.
- > Programmable lifting speeds
- > PDA handset or standard wired hand control
- > Downloadable/printable compliance software
- > Down loadable/printable service software
- > Integrated lift color display
- > Blue ambient light – on top of hoist
- > Multiple spreader bar options or -
- > Electronic Patient Positioner
- > Approved for use in wet and dry rooms
- > Horizontal drive motor
- > Soft start/stop function
- > Quick Release hanger Bar
- > Permanent Charged rail

UNILIFT/PRO

8 function Medical PDA

Used as a hand control, read data, remote down load.

- > 5" high resolution touch screen
- > Anti-bacterial material
- > Readable and downloadable compliance data
- > Readable and downloadable service data
- > Wet hands operating is available
- > Charging station
- > Drop resistance
- > Built in user instructions with illustrations via PDA




Hand Control



Innovative rail solutions

Supporting architects and healthcare facilities by providing the ideal track solution, meeting the clinical need, while maintaining a contemporary aesthetic look. With 32 years of collaborating with healthcare professionals, architects, and patients you can be confident in choosing barrier free lifts.

Adjustable switches guide the lift into the different guiderails. Additionally, in traverse systems, gates ensure expansion of the track layout.



A: Carousel

B: Track Switch

C: Gate, Linear track to transverse rail

D: Transverse System

