



Freeslides

Slide Tubes & Slide Sheets

Single Patient Use

Freeslides can be used for a multitude of patient handling tasks. The anti-friction tube and sheets reduce skin shear making transfers safer and easier for staff and patients. Patient comfort is increased while staff muscular skeletal injuries are decreased. All Freeslides are designed to be wiped down between each use. They are made from 100% polyester with an Anti-Static polyurethane coating.

Features:

- **Wall Dispenser:** Readily available
- **Heat Seal Edges:** No seams to create friction. New!
- **Versatile:** Can be used for more than 25 SPH tasks
- **Single Patient Use:** Reduces the risk of cross-contamination
- **Cost Effective :** Safe transfers in any locations, on a budget
- **No Laundering:** No loss, no logistical challenges with availability.
- **Always There:** Travels from unit to unit throughout patients stay
- **OSHA & ACON:** Recognized Tubes and Sheets for SPH
- **Recyclable:** Recyclable material



Slide Tubes: : One Sheet Tied Into A Continuous Loop

Sizes: 58" X 48" and 58" X 29"



Slide Sheets: Flat Sheet

Sizes: 58" X 48" and 58" X 29"

